Description

Evetric is a light-curing, radiopaque nano-hybrid composite for restorative treatment. Evetric cures with light in the wavelength range of 400–500 nm (blue light).

Composition

Evetric consists of dimethacrylates (19-20 wt.%). The fillers contain barium glass, ytterbium trifluoride, mixed oxide and copolymers (80-81 wt.%). Additives, catalysts, stabilizers and pigments are additional contents (< 1 wt.%). The total content of inorganic fillers is 55–57 vol.%. The particle size of inorganic fillers is between 40 nm and 3000 nm.

Indications

- Class I-V restorations
- Restoration of deciduous teeth
- Extended fissure sealing
- Direct veneers
- Splinting of mobile teeth
- Blocking out of undercuts
- Repair of composite and ceramic veneers

Contraindications

- If a dry working field cannot be established or the stipulated application technique cannot be applied.
- If the patient is known to be allergic to any of the ingredients.

Side effects

- In rare cases, components of Evetric may lead to sensitization. The product should not be used in such cases.
- To avoid possible irritation of the pulp, areas close to the pulp should be protected with a suitable pulp/dentin protector.

Interactions

- Materials containing eugenol or clove oil may inhibit the polymerization of Evetric. Consequently, the application of such materials together with the product must be avoided.
- Discolouration may occur in combination with cationic mouthwashes, plaque disclosing agents and chlorhexidine.

Application**Direct Restorations****1. Shade determination**

Clean the teeth prior to shade determination. The shade is selected with the tooth still moist.

2. Isolation

Appropriate isolation, best with a rubber dam (e.g. OptraDam® Plus), is required.

3. Cavity preparation

Prepare the cavity according to the requirements of the adhesive technique.

Subsequently, remove all residue in the cavity with water spray and dry with water- and oil-free air.

4. Pulp protection / Base

Do not apply a base material when using an enamel/dentin bonding agent. Only cover very deep areas close to the pulp with a calcium hydroxide material (e.g. ApexCal®) and subsequently use a pressure-resistant cement (e.g. a glass ionomer cement, such as Vivaglass® Liner). Do not cover other cavity walls, since they can be used to support the bond with an enamel/dentin adhesive.

5. Matrix / Interdental wedge

Place a matrix and an interdental wedge, if required.

6. Conditioning / Application of the bonding agent

Condition and apply the bonding agent according to the Instruction for Use of the product in use. We recommend using Syntac® (with phosphoric acid etching) or ExcITE® F (with phosphoric acid etching) or the self-etching adhesive AdheSE® One F.

7. Application of Evetric

Apply Evetric in layers of max. 2mm or 1.5mm (Dentin shade) and contour/adapt the material to the cavity walls using a suitable instrument (e.g. OptraSculpt®). Use a special contact point instrument (e.g. OptraContact) for large cavities where the contacts are difficult to create. Polymerize each layer individually according to the table below. Hold the light emission window as closely as possible to the surface of the restorative material.

A flowable composite can be used as a thin initial layer. This optional layer has to be cured separately (please refer to the respective Instruction for Use).

8. Finishing / Checking the occlusion / Polishing

Finish the restoration with suitable finishers or fine diamonds. Subsequently, check the occlusion. Polish with silicone polishers (e.g. OptraPol® Next Generation).

Unit Program	LEdition	bluephase C8	bluephase	bluephase 20i
Turbo	—	—	—	5 s
High Power	20 s	15 s	10 s	10 s
Soft Start	—	20 s	15 s	15 s

Light intensity	Exposure Time
≥ 500 mW/cm ²	20 s
≥ 1,000 mW/cm ²	10 s

Special note

1. If fresh Evetric is applied to polymerized material, it must first be roughened and wetted with Heliobond.
2. Evetric should have room temperature when applied. Cool temperatures render the material difficult to extrude.
3. Do not use oxidizing disinfectants to disinfect syringes.

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– Light-curing resin-based dental restorative material

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Warning

- Avoid contact of unpolymerized material with the skin, mucous membrane and eyes.
- Unpolymerized material may have a slight irritating effect and may lead to a sensitization against methacrylates.
- Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

Storage and hygiene

- Storage temperature: 2-28 °C (36-82 °F)
- Close syringes immediately after usage. Exposure to light causes premature polymerization.
- Expiry date: see note on syringes and packaging.
- Do not use any of these products after the indicated date of expiration.

Keep out of the reach of children.

For use in dentistry only.

This material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damage resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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Ivoclar Vivadent Marketing Ltd.
Ivoclar Vivadent AB
Ivoclar Vivadent Limited
Ivoclar Vivadent, Inc.